From 50mT to 7T: *Pulseq* Live Across Hardware Platforms and Field Strengths

Maxim Zaitsev¹, Thomas O'Reilly², Benjamin Menkuec³, Marcus Prier⁴, David Schote⁴, Jason Stockmann^{5,6}, Lincoln Craven-Brightman^{5,6}, Melissa Haskell^{7,8}, Jon-Fredrik Nielsen⁸

- 1. High Field Magnetic Resonance Center, Center for Medical Physics and Biomedical Engineering, Medical University of Vienna, Vienna, Austria
- 2. Department of Radiology, Leiden University Medical Center, Leiden, Netherlands
- 3. University of Applied Sciences and Arts Dortmund, Dortmund, Germany
- 4. Otto von Guericke Universität Magdeburg, Magdeburg, Germany
- 5. Department of Radiology, Massachusetts General Hospital, Charlestown, MA, USA
- 6. Harvard Medical School, Boston, MA, USA
- 7. Electrical Engineering and Computer Science, University of Michigan, Ann Arbor, MI, USA
- 8. fMRI laboratory and Biomedical Engineering, University of Michigan, Ann Arbor, MI, USA.









Outline


- *Pulseq* sequence definition and programming language overview
- Using *Pulseq* on different types of scanners:
 - Siemens, GE, and tabletop scanners using OCRA
- Live demo of *Pulseq* across sites
 - 1. Basic pulse-acquire scans: Free Induction Decay (FID), Spin-Echo (SE)
 - 2. Cartesian spin-echo imaging
 - 3. Radial imaging (spin-echo and gradient echo)
 - 4. Advanced topics (gradient delays and correction, etc...)





Pulseq sequence programming

- Cross-platform pulse programming framework for MRI
- Low level: Pulseq file
- High level: MATLAB¹ or Python² toolboxes



Source: Layton, et al., "Pulseq: A rapid and hardware-independent pulse sequence prototyping framework", MRM 2017

- 1. http://pulseq.github.io
- 2. http://github.com/imr-framework/pypulseq






Pulseq file

• Explicit (low level) specification of the pulse sequence

Think of an MP3 file (or more precisely lossless FLAC)

No loops, no parameters, no dependencies, no fuss!

- Text file (human-readable)
 - Simple hierarchy (RF pulses, gradients, shapes)
 - Event table keeps it together
 - See http://pulseq.github.io/specification.pdf
 for more details








Pulse sequence definition

Concatenation of non-overlapping blocks



- Block 1: gradient and RF
- Block 2: only gradients
- Block 3: gradient and ADC
- Block 4: only gradients
- Block 5: gradient and RF ...

Gradients do not have to start or end at 0 at the block boundaries (advanced feature available on some hardware platforms)

High-level programming environments

- Matlab Pulseq toolbox
- Python *pypulseq* toolbox





- Further options
 - Pulseq-GPI a graphical sequence programming environment
 - TOPPE is primarily targeted at GE but can import and export *pulseq* files
 - GammaStar can export *pulseq* files
 - JEMRIS Bloch simulator can export *pulseq* files
 - CoreMRI Bloch simulator can export *pulseq* files

• ...





Matlab Pulseq workflow

- Define the system properties
- Define high-level parameters (convenience)
- Define pulses used in the sequence
- Calculate the delays and reordering tables
- Loop and define sequence blocks
 - Duration of each block is defined by the duration of the longest event
- Copy 'gre.seq' to the scanner and run it!
- Screenshot shows an entire runnable gradient echo sequence code (similar to Siemens' example miniFlash)

```
system = mr.opts('MaxGrad',30,'GradUnit','mT/m',...
    'MaxSlew',170, 'SlewUnit', 'T/m/s');
seq=mr.Sequence(system);
fov = 220e-3; Nx=64; Ny=64; TE = 10e-3; TR = 20e-3;
[rf, gz] = mr.makeSincPulse(15*pi/180, system, 'Duration', 4e-3,...
    'SliceThickness', 5e-3, 'apodization', 0.5, 'timeBwProduct', 4);
qx = mr.makeTrapezoid('x', system, 'FlatArea', Nx/fov, 'FlatTime', 6.4e-3);
adc = mr.makeAdc(Nx,'Duration',qx.flatTime,'Delay',qx.riseTime);
gxPre = mr.makeTrapezoid('x', system, 'Area', -gx.area/2, 'Duration', 2e-3);
gzReph = mr.makeTrapezoid('z',system,'Area',-gz.area/2,'Duration',2e-3);
phaseAreas = ((0:Ny-1)-Ny/2)*1/fov;
delayTE = TE - mr.calcDuration(gxPre) - mr.calcDuration(rf)/2 ...
    - mr.calcDuration(gx)/2;
delayTR = TR - mr.calcDuration(gxPre) - mr.calcDuration(rf) ...
    - mr.calcDuration(gx) - delayTE;
delay1 = mr.makeDelay(delayTE);
delay2 = mr.makeDelay(delayTR);
for i=1:Nv
    seq.addBlock(rf,gz);
    gyPre = mr.makeTrapezoid('y',system,'Area',phaseAreas(i),...
                              'Duration', 2e-3);
    seq.addBlock(gxPre,gyPre,gzReph);
    seq.addBlock(delay1);
    seq.addBlock(qx,adc);
    seq.addBlock(delay2)
seq.write('gre.seg')
```






Development and debugging

- Visualize sequences with seq.plot()
- Verify raster alignment and system delays with seq.checkTiming()
- Check echo conditions and visualize k-spaces with seq.calculateKspace() and calculateKspacePP()
- Check details of gradient waveforms with seq.gradient_waveforms()
- Evaluate further details (TE,TR, slew rate, etc...) with seq.testReport()
- Listen to seq.sound()





Plotting the sequence timing




• Further plot options: 'timeRange', 'TimeDisp', 'showBlocks', 'Label'


Pulseq unit conventions

- Default to SI units where possible
- Default gradient unit is Hz/m
 - Other units possible
- K-space unit is 1/m
- FOV calculations are trivial

$$k_{max} = 1/\Delta x \quad \Delta k = 1/FOV$$

- No γ required!
- No 2π required!
- K-space analysis is trivial
 - Spatial resolution of the sequence in Fig. 2 is 1m/125=4mm







- PyPulseq is a close replica of the original Pulseq toolbox that does not require a MATLAB license
- Runs in many of established Python environments,
 e.g. in Jupyter notebook (http://jupyter.org/)
- PyPulseq may now run in your browser via the Google Colaboratory
- To learn more about *PyPulseq* check out ISMRM 2021 e-poster #3760 on Thursday or visit http://github.com/imr-framework/pypulseq



```
pypulseq_test.ipyn ×
 C û □ □ https://colab.research.google.com/dri ··· □ ☆
  ♠ pypulseq test.ipynb ☆
  File Edit View Insert Runtime Tools Help Last edited on A
     # (re)initialize the sequence object
      seg = Seguence(system)
      # Prepare RF offsets. This is required for multi-slice acquisition
      delta_z = n_slices * slice_gap
      z = np.linspace((-delta_z / 2), (delta_z / 2), n_slices) + rf_offset
      for k in range(nsa): # Averages
        for j in range(n slices): # Slices
         freq offset = gz90.amplitude * z[j]
          rf90.freq_offset = freq_offset
         rf180.freq offset = freq offset
          for i in range(Nv): # Phase encodes
            seq.add block(rf90, gz90)
            gv pre = make trapezoid(channel='v', system=system.
                                  area=phase_areas[-i -1], duration=2e-3)
            seg.add block(delav1)
            seg.add block(gz spoil)
            seq.add block(delay2)
            seq.add block(qx, adc)
            gy_pre = make_trapezoid(channel='y', system=system,
                                  area=-phase_areas[-j -1], duration=2e-3)
 [ ] seq.plot(time range=(0, 0.12))
                 0.02 0.04 0.06
                                     0.08
                 0.02 0.04 0.06 0.08 0.10 0.12
                 0.02 0.04 0.06 0.08 0.10 0.12
[ ] seq.check_timing()
```



Outline


- Pulseq sequence definition and programming language overview
- Using *Pulseq* on different types of scanners:
 - Siemens, GE, and tabletop scanners using OCRA
- Live demo of *Pulseq* across sites
 - 1. Basic pulse-acquire scans: Free Induction Decay (FID), Spin-Echo (SE)
 - 2. Cartesian spin-eco imaging
 - 3. Radial imaging (spin-echo and gradient echo)
 - 4. Advanced topics (gradient delays and correction, etc...)





Pulseq on Siemens scanners

- Optional step: connect an additional PC to the scanner
- Save the .seq file on the scanner as external.seq
- Run the interpreter_sequence on the scanner
- Optional step: stream raw data to your PC with NIH_DataCatcher
- or export raw data manually









IDEA simulation with *Pulseq*

Pulseq interpreter sequence can also be used with the Siemens' IDEA

- 1. Save your .seq file as %CustomerSeq%/Pulseq/external.seq
- In the IDEA command run sim







To acquire the Pulseq interpreter sequence contact pulseq.mr@uniklinik-freiburg.de





Pulseq on GE scanners

Step 1: Download the 'PulseqGEq' and 'TOPPE' Matlab toolboxes

```
$ git clone git@github.com:toppeMRI/PulseGEq.git
$ git clone git@github.com:toppeMRI/toppe.git
```

Step 2: Set system hardware parameters (max gradient, slew, etc)


```
>> sys = toppe.systemspecs('maxGrad',5,'gradUnit','G/cm');
```

Step 3: Convert the .seq file to the TOPPE file format


```
>> pulsegeq.seq2ge('se radial.seq','system',sys);
```

Step 4: Preview the sequence in movie/loop mode

```
>> toppe.playseq(3);
```



Step 5: Copy 'toppeScanFiles.tar' to the scanner, untar, and run the TOPPE interpreter sequence








Pulseq on tabletop scanners using OCRA/FLOCRA

OCRA/FLOCRA: Open source consoles compatible with tabletop and custom scanners


- Runs on Red Pitaya 122 (FLOCRA/OCRA) and 125 (OCRA) boards
- Open source Python code runs MR experiments using **Pulseq** .seq files



0.48T Tabletop Dortmund, Germany



0.35T Tabletop Boston, MA, USA



50mT Halbach Leiden, Netherlands



Red Pitaya 122 board

FLOCRA: \$ git clone https://github.com/lcbMGH/flocra-pulseq.git

OCRA: \$ git clone --recurse-submodules https://github.com/OpenMRI/ocra-pack.git






Pulseq on tabletop scanners using OCRA/FLOCRA

RP-122/125 Board Hardware:

- Two RF inputs, two RF outputs, capable of 300 kHz to 100 MHz signals.
- 16 digital pins for extensions used for gradient SPI bus and TR-switch control
- 125 (RP-125), 122.88 (RP122) MS/s sample rate
- Compatible with OCRA1 and GPA-FHDO gradient DAC boards
- Ethernet connection for control



FLOCRA: \$ git clone https://github.com/lcbMGH/flocra-pulseq.git

OCRA: \$ git clone --recurse-submodules https://github.com/OpenMRI/ocra-pack.git

Outline







- *Pulseq* sequence definition and programming language overview
- Using *Pulseq* on different types of scanners:
 - Siemens, GE, and tabletop scanners using OCRA
- Live demo of *Pulseq* across sites
 - 1. Basic pulse-acquire scans: Free Induction Decay (FID), Spin-Echo (SE)
 - 2. Cartesian spin-eco imaging
 - 3. Radial imaging (spin-echo and gradient echo)
 - 4. Advanced topics (gradient delays and correction, etc...) are only covered in the live portion of the tutorial



- Example 1: free induction decay sequence
- After line 25
 the sequence is ready
 (in Matlab's memory)
- Further lines: display, sanity checks, export to disk



 Example 2: FID sequence with flip angle variation



```
/home/zaitsev/pulseq_home/pulseq_public/matlab/ismrmTutorial/seq/seq02_Fid_multipleFAs.m -
                                                                                        EDITOR
                 PUBLISH
                 Find Files
                               ♦
                                                                         Run Section
                 Compare
                               Go To ▼
                                             Breakpoints
                                                                 Run and
                                                                         Advance
                A Print -
                               Q Find ▼
                                                                                        Time
                                                                 Advance
                               NAVIGATE
                                             BREAKPOINTS
        system = mr.opts('rfRingdownTime', 20e-6, 'rfDeadTime', 100e-6, ...
 2
                         'adcDeadTime', 20e-6);
 3
4 -
       seq=mr.Sequence(system);
                                              % Create a new sequence object
       Nx=4096;
       Nrep=1;
       adcDur=51.2e-3:
       rfDur=1000e-6; % increased to 1ms
       TR=2000e-3;
                     % may increase to ~5s avoid T1 saturation
11 -
       flip angles=18;
12
13
       % Define delays and ADC events
14 -
       adc = mr.makeAdc(Nx, 'Duration', adcDur, 'system', system, 'delay', TE-rfDur/2-system.rfRingdown
15
16 -
       delavTR=TR-mr.calcDuration(rf):
17 -
       assert(delavTR>=0);
18
19
       % Loop over repetitions and define sequence blocks
20
      □ for f=1:flip angles
21
            % Create non-selective pulse
22 -
            rf = mr.makeBlockPulse(pi/flip angles*f,'Duration',rfDur, 'system', system');
            % Loop over repetitions and define sequence blocks
23
24 -
            for i=1:Nrep
25 -
                seq.addBlock(rf);
26 -
                seq.addBlock(adc,mr.makeDelay(delayTR));
27 -
28 -
       end
29
30
       % show the entire sequence
31 -
       seq.plot();
32
33
       % check whether the timing of the sequence is compatible with the scanner
34 -
       [ok, error report]=seq.checkTiming;
35
36
       if (ok)
37 -
            fprintf('Timing check passed successfully\n');
38 -
39 -
            fprintf('Timing check failed! Error listing follows:\n');
40 -
            fprintf([error report{:}]);
            fprintf('\n');
41 -
42 -
       end
43
44 -
       seq.write('fid.seq')
                                   % Write to pulseq file
45
       %seq.install('siemens');
                                   % copy to scanner
                                            UTF-8
                                                             script
                                                                                                        Col 1
                                                                                                 Ln 3
```

- Example 3: naive spin-echo
- The sequence is ready after line 34
- We (ab)use k-space calculation routines for their ability to deliver timestamps

two RF pulses

check TE calculation

```
/home/zaitsev/pulseq_home/pulseq_public/matlab/ismrmTutorial/seq/seq03_SE.m
                Tind Files
                                           Insert 🗐 fx 📆
                                                                                Run and 🖳 Advance
       system = mr.opts('rfRingdownTime', 20e-6, 'rfDeadTime', 100e-6,
                         'adcDeadTime', 20e-6);
       seq=mr.Sequence(system);
                                             % Create a new sequence object
       Nrep=1:
       adcDur=51.2e-3;
       rfDur=1000e-6;
       TR=100e-3:
       TE=54e-3;
11
12
       % Create non-selective excitation and refocusing pulses
13 -
       rf ex = mr.makeBlockPulse(pi/2,'Duration',rfDur, 'system', system);
       rf_ref = mr.makeBlockPulse(pi,'Duration',rfDur, 'system', system', 'PhaseOffset',pi/2, 'use', 'refocusing');
14 -
15
16
       % Define delays and ADC events
17
       delayTE1=TE/2-mr.calcDuration(rf_ex)/2-mr.calcDuration(rf_ref)/2;
       delayTE2=TE/2-mr.calcDuration(rf ref)+rf ref.delay+mr.calcRfCenter(rf ref)-adcDur/2; % this is not perfect, bu
19
20
       adc = mr.makeAdc(Nx, 'Duration', adcDur, 'system', system, 'delay', delayTE2);
21
       delayTR=TR-mr.calcDuration(rf ex)-delayTE1-mr.calcDuration(rf ref);
23
       assert(delayTE1>=0);
       assert(delayTE2>=0);
       assert(delavTR>=0):
       % Loop over repetitions and define sequence blocks
           seq.addBlock(rf ex);
           seq.addBlock(mr.makeDelay(delayTE1));
           seq.addBlock(rf ref):
           seq.addBlock(adc,mr.makeDelay(delayTR));
34
       % show the entire sequence showing the block structure
       seq.plot('showBlocks',1,'timeDisp','ms');
       % check whether the timing of the sequence is compatible with the scanner
       [ok, error_report]=seq.checkTiming;
42
           fprintf('Timing check passed successfully\n');
           fprintf('Timing check failed! Error listing follows:\n');
           fprintf([error_report{:}]);
47
           fprintf('\n');
48
49
50
       seq.write('se.seq')
                                 % Write to pulseq file
       %seq.install('siemens');
52
53
       % calculate k-space but only use it to check the TE calculation
54 -
       [ktraj adc, t adc, ktraj, t ktraj, t excitation, t refocusing] = seg.calculateKspacePP();
55
       assert(abs(t_refocusing-t_excitation-TE/2)<le-6); % check that the refocusing happens at the 1/2 of TE
57 -
       assert(abs(t adc(Nx/2)-t excitation-TE)<adc.dwell); % check that the echo happens as close as possible to the
                                                    UTF-8
                                                                                                        Ln 3
```

- Example 4: spin-echo with spoilers around the refocusing pulse
- Create spoiling gradient
- Use it in the sequence in two different blocks



```
/home/zaitsev/pulseq_home/pulseq_public/matlab/ismrmTutorial/seq/seq04_SE_withSpolers.m
                Tind Files
                              BREAKPOINTS
                                                  rfDeadTime', 100e-6,
       seq=mr.Sequence(system);
                                             % Create a new sequence object
       adcDur=51.2e-3;
       rfDur=1000e-6;
       TR=250e-3;
       spA=1000; % spoiler area in 1/m (=Hz/m*s)
       % todo: change spoiler area, remove one of spoilers and observe the signal
       % Create non-selective excitation and refocusing pulses
       rf ex = mr.makeBlockPulse(pi/2, 'Duration', rfDur, 'system', system);
       rf_ref = mr.makeBlockPulse(pi,'Duration',rfDur, 'system', system, 'PhaseOffset' ,pi/2 , 'use', 'refocusing');
       % calculate spoiler gradient, let's put it on X axis for now
       g_sp=mr.makeTrapezoid('x','Area',spA,'system',system);
       rf_ref.delay=max(mr.calcDuration(g_sp),rf_ref.delay);
       % Define delays and ADC events
       delayTE1=TE/2-(mr.calcDuration(rf ex)-mr.calcRfCenter(rf ex)-rf ex.delay)-rf ref.delay-mr.calcRfCenter(rf ref);
       delayTE2=TE/2-mr.calcDuration(rf_ref)+rf_ref.delay+mr.calcRfCenter(rf_ref)-adcDur/2; % this is not perfect, but
       assert(delayTE2>mr.calcDuration(g sp));
       adc = mr.makeAdc(Nx,'Duration',adcDur, 'system', system, 'delay', delayTE2);
       delayTR=TR-mr.calcDuration(rf_ex)-delayTE1-mr.calcDuration(rf_ref);
       assert(delavTE1>=0);
       assert(delayTE2>=0);
       assert(delayTR>=0);
       % Loop over repetitions and define sequence blocks
36 -
      37 -
           seq.addBlock(rf_ex);
38 -
           seq.addBlock(mr.makeDelay(delayTE1));
39
           seq.addBlock(rf_ref,g_sp);
40 -
           seq.addBlock(adc,g sp,mr.makeDelay(delayTR));
41 -
       seq.plot(); % show the entire sequence
       % check whether the timing of the sequence is compatible with the scanner
       [ok, error_report]=seq.checkTiming;
49 -
           fprintf('Timing check passed successfully\n');
50
51 -
           fprintf('Timing check failed! Error listing follows:\n');
52
           fprintf([error_report{:}]);
53 -
           fprintf('\n');
54
       seq.write('se.seq')
                                 % Write to pulseq file
       %seq.install('siemens'); % copy to scanner
59
       % calculate k-space but only use it to check timing
60 -
       [ktraj adc, t adc, ktraj, t ktraj, t excitation, t refocusing] = seq.calculateKspacePP();
       assert(abs(t refocusing-t excitation-TE/2)<le-6); % check that the refocusing happens at the 1/2 of TE
       assert(abs(t adc(Nx/2)-t excitation-TE)<adc.dwell); % check that the echo happens as close as possible to the mi
                                                                                                         Ln 43 Col 39
```


• Example 5: Cartesian spin-echo imaging sequence

dummy

calculate

delays

- Use shorter ADC window compared to the spectroscopy sequence
- More gradients to create
- More complex delay calculation
- The sequence is ready at line 67
- See timing checks in the live portion








Pulseq - seq05 contd.

- Example 5:
 Cartesian spin-echo imaging sequence
- Very useful option in the Pulseq toolbox is the k-space calculation
- See how we travel in the k-space to the bottom right in the last TR to be "teleported" to the upper left corner (see also k(t) plots in the live session)



- Example 6:2D radial spin-echo imaging sequence
- Remove phase encoding
- Rotate read prephaser and readout gradients
- The sequence is actually simpler
- The sequence is ready at line 63
- See timing checks in the live portion


```
/home/zaitsev/pulseq_home/pulseq_public/matlab/ismrmTutorial/seq/seq06_RadialSE.m
                             4
                                          Insert 🛃 fx 👍
       seg=mr.Seguence(system);
                                             % Create a new sequence object
       adcDur=2.56e-3:
       rfDur1=3e-3:
       rfDur2=le-3;
       TR=150e-3;
       TE=54e-3;
       spA=1000; % spoiler area in 1/m (=Hz/m*s) % MZ: need 5000 for my oil phantom
       sliceThickness=3e-3;
       fov=250e-3: Nx=256:
                                       % Define FOV and resolution
                                       % number of radial spokes
       Ndummy=10;
                                       % number of dummy scans
       delta=pi/Nr;
                                       % angular increment; try golden angle pi*(3-5^0.5) or 0.5 of it
                                       % the orientation of the first projection
       % Create 90 degree slice selection pulse and gradient and even the
       % refocusing gradient; we will not use it however but will subtract its
       % area from the first spoiler
       [rf ex, gs, gsr] = mr.makeSincPulse(pi/2,system,'Duration',3e-3,...
             SliceThickness',sliceThickness,'apodization',0.4,'timeBwProduct',4);
       % Create non-selective refocusing pulse
       rf ref = mr.makeBlockPulse(pi, 'Duration', rfDur2, 'system', system', 'PhaseOffset', pi/2, 'use', 'refocusing');
       % calculate spoiler gradient
       g_spl=mr.makeTrapezoid('z','Area',spA+gsr.area,'system',system);
       rf_ref.delay=max(mr.calcDuration(g_sp1),rf_ref.delay);
       g sp2=mr.makeTrapezoid('z','Area',spA,'system',system);
       % Define delays and ADC events
       deltak=1/fov;
       gr = mr.makeTrapezoid('x',system,'FlatArea',Nx*deltak,'FlatTime',adcDur);
       adc = mr.makeAdc(Nx,system,'Duration',adcDur,'delay',gr.riseTime);
       grPredur = 5e-3; % use a fixed time to make this gradient visible on the plot
       grPre = mr.makeTrapezoid('x',system,'Area',gr.area/2+deltak/2,'Duration',grPredur);
       delayTE1=TE/2-(mr.calcDuration(gs)-mr.calcRfCenter(rf ex)-rf ex.delay)-rf ref.delay-mr.calcRfCenter(rf ref)-grF
       delayTE2=TE/2-mr.calcDuration(rf_ref)+rf_ref.delay+mr.calcRfCenter(rf_ref)-adc.delay-adcDur/2; % this is not pe_
       delayTR=TR-mr.calcDuration(gs)-grPredur-delayTE1-mr.calcDuration(rf_ref)-delayTE2-mr.calcDuration(gr);
       assert(delayTE1>=0);
        assert(delayTE2>mr.calcDuration(g sp2));
       assert(delayTR>=0);
       % Loop over repetitions and define sequence blocks
        for i=(1-Ndummy):Nr
           seq.addBlock(rf_ex,gs);
                                                                      Rotate #1
           seq.addBlock(mr.rotate('z',delta*(i-1),grPre));
           seq.addBlock(mr.makeDelay(delayTE1));
           seq.addBlock(rf_ref,g_sp1);
           seq.addBlock(g_sp2, mr.makeDelay(delayTE2));
58 -
59 -
60 -
               seq.addBlock(mr.rotate('z',delta*(i-1),adc,gr));
                                                                             Rotate #2
               seq.addBlock(mr.rotate('z',delta*(i-1),gr));
61 -
62 -
63 -
64
           seq.addBlock(mr.makeDelay(delayTR));
       seq.plot();
                                                                                                           Ln 3
```






Pulseq – seq06 contd.

Time to explore sequence timing and the trajectory





- Example 7: 2D radial gradient echo
- Add RF spoiling
 - Introduce counters
 - Update phase of the RF pulse
 - Remember to update ADC phase
- Add gradient spoiler at the end of the readout gradient
- The sequence is ready at line 63
- See timing checks in the live portion

gradient spoiling

```
/home/zaitsev/pulseq_home/pulseq_public/matlab/ismrmTutorial/seq/seq07_RadialGradientEcho.m
         Find Files
                      sea=mr.Seauence():
                                % Create a new sequence object
fov=250e-3: Nx=256:
                                % Define FOV and resolution
alpha=30;
                                % flip angle
sliceThickness=3e-3;
                                % TE; give a vector here to have multiple TEs (e.g. for field mapping)
TR=20e-3:
                                % only a single value for now
Nr=256;
                                % number of radial spokes
                                % number of dummy scans
Ndummy=10;
                                % angular increment; try golden
                                                                      pi*(3-5^0.5) or 0.5 of it
delta=pi/Nr;
% more in-depth parameters
                                                                   RF spoiling
rfSpoilingInc=117;
                                % RF spoiling increment
% set system limits
sys = mr.opts('MaxGrad', 28, 'GradUnit', 'mT/m', ...
    'MaxSlew', 80, 'SlewUnit', 'T/m/s', 'rfRingdownTime', 20e-6, .
    'rfDeadTime', 100e-6, 'adcDeadTime', 10e-6);
% Create alpha-degree slice selection pulse and gradient
[rf, gz] = mr.makeSincPulse(alpha*pi/180,'Duration',4e-3,...
     SliceThickness', sliceThickness, 'apodization', 0.5, 'timeBwProduct', 4, 'system', sys);
% Define other gradients and ADC events
gx = mr.makeTrapezoid('x', 'FlatArea', Nx*deltak, 'FlatTime', 6.4e-3, 'system', sys);
adc = mr.makeAdc(Nx,'Duration',qx.flatTime,'Delay',qx.riseTime,'system',sys);
gxPre = mr.makeTrapezoid('x','Area',-gx.area/2-deltak/2,'Duration',2e-3,'system',sys); % we need this "d
gzReph = mr.makeTrapezoid('z','Area',-gz.area/2,'Duration',2e-3,'system',sys);
% gradient spoiling
gxSpoil=mr.makeTrapezoid('x','Area',0.5*Nx*deltak,'system',sys);
qzSpoil=mr.makeTrapezoid('z','Area',4/sliceThickness,'system',sys);
% Calculate timing
delayTE=ceil((TE - mr.calcDuration(qxPre) - qz.fallTime - qz.flatTime/2 ...
      mr.calcDuration(gx)/2)/seq.gradRasterTime)*seq.gradRasterTime;
delayTR=ceil((TR - mr.calcDuration(gxPre) - mr.calcDuration(gz) ...

    mr.calcDuration(gx) - delayTE)/seq.gradRasterTime)*seq.gradRasterTime;


assert(all(delayTR>=mr.calcDuration(gxSpoil,gzSpoil)));
rf inc=0;
                                                                           RF
for i=(-Ndummy):Nr
    for c=1:length(TE)
        rf.phaseOffset=rf phase/180*pi;
                                                                     spoiling
        adc.phaseOffset=rf_phase/180*pi;
        rf inc=mod(rf inc+rfSpoilingInc, 360.0);
        rf phase=mod(rf phase+rf inc, 360.0);
        seq.addBlock(rf,gz);
        phi=delta*(i-1);
        seq.addBlock(mr.rotate('z',phi,gxPre,gzReph));
        seq.addBlock(mr.makeDelay(delayTE(c)));
            seq.addBlock(mr.rotate('z',phi,gx,adc));
            seq.addBlock(mr.rotate('z',phi,gx));
        seq.addBlock(mr.rotate('z',phi,gxSpoil,gzSpoil,mr.makeDelay(delayTR)));
seq.plot();
                                                                                            Ln 10 Col 1
```






Pulseq - seq07 contd.

Time to explore sequence timing and the trajectory
 (e.g. RF phase and the gradient spoiling effects on the trajectory)






Pulseq – recon scripts

file number selector

- All recon scripts are similar
- General logic
 - We expect data to be stored in a directory of choice (line 7 in this example points to data)
 - Raw data MUST always be accompanied with a Pulseq file with the identical name
 - By default all scripts load the most recent data
 - To load the second last data set replace 'end-0' with 'end-1' (line 14)
 - Ideally we expect raw data in a 3D array with dimensions: samples x channels x readouts
 - All data counters, data sorting and similar parameters is calculated from the Pulseq file
- Here is the complete FID/SE example



Pulseq – 2D FFT


k-space trajectory

- The beginning of the script is similar to the FID/SE example
- Screen shot shows automatic counter-free data sorting part
- k-space trajectory automatically updates data reordering for (almost) any 2D Cartesian sequence
- Following (not shown) part is a trivial 2D FFT and a sum-of-square combination for multichannel data

```
/home/zaitsev/pulseq_home/pulseq_public/matlab/ismrmTutorial/recon/recon02_2DFFT.m
                                          Insert 💂 fx 👍 ▼
        [adc len,channels,readouts]=size(data unsorted);
       % the incoming data order is [kx coils acquisitions]
       data coils last = permute(data unsorted, [1, 3, 2]);
       % Plot and analyze the trajectory data (ktraj_adc)
       [ktraj_adc, ktraj, t_excitation, t_refocusing, t_adc] = seq.calculateKspace();
50 -
       figure; plot(ktraj(1,:),ktraj(2,:),'b',...
                    ktraj_adc(1,:),ktraj_adc(2,:),'r.'); % a 2D plot
       axis('equal'); title('2D k-space trajectory');
53
       % try to detect the data ordering
       k_extent=max(abs(ktraj_adc),[],2);
       k scale=max(k extent);
        k threshold=k scale/5000;
                                                                              remove slice
       % detect unused dimensions and delete them
       if any(k extent<k threshold)
           ktraj adc(k extent<k threshold,:)=[]; % delete rows
                                                                               dimension
           k extent(k extent<k threshold)=[];
       % detect dK, k-space reordering and repetitions (or slices, etc)
       kt sorted=sort(ktraj adc,2);
       dk_all=kt_sorted(:,2:end)-kt_sorted(:,1:(end-1));
        dk_all(dk_all<k_threshold)=NaN;
                                                                                         calculate
        dk min=min(dk all,[],2);
        dk max=max(dk all,[],2);
        dk_all(dk_all-dk_min(:,ones(1,size(dk_all,2)))>k_threshold)=NaN;
                                                                                       reordering
        dk_all_cnt=sum(isfinite(dk_all),2);
        dk_all(~isfinite(dk_all))=0;
        dk=sum(dk_all,2)./dk_all_cnt;
        [~,k0 ind]=min(sum(ktraj adc.^2,1));
        kindex=round((ktraj_ado-ktraj_ado(:,k0_ind*ones(1,size(ktraj_ado,2))))./dk(:,ones ,size(ktraj_ado,2))));
        kindex min=min(kindex,[],2);
       kindex mat=kindex-kindex min(:,ones(1,size(ktraj adc,2)))+1;
       kindex end=max(kindex mat,[],2);
        sampler=zeros(kindex end');
        repeat=zeros(1,size(ktraj_adc,2));
        for i=1:size(kindex_mat,2)
           if (size(kindex mat,1)==3)
               ind=sub2ind(kindex end,kindex mat(1,i),kindex mat(2,i),kindex mat(3,i));
               ind=sub2ind(kindex_end,kindex_mat(1,i),kindex_mat(2,i));
           repeat(i)=sampler(ind);
           sampler(ind)=repeat(i)+1:
       if (max(repeat(:))>0)
           kindex=[kindex;(repeat+1)];
           kindex_mat=[kindex_mat;(repeat+1)];
           kindex_end=max(kindex_mat,[],2);
       %figure; plot(kindex(1,:),kindex(2,:),'.-');
        % sort the k-space data into the data matrix
       data_coils_last=reshape(data_coils_last, [adc_len*readouts, channels]);
                                                                                    sort the data
        data=zeros([kindex_end' channels]);
       if (size(kindex,1)==3)
           for i=1:size(kindex.2)
               data(kindex mat(1,i),kindex mat(2,i),kindex mat(3,i),:)=data coils
           for i=1:size(kindex,2)
               data(kindex_mat(1,i),kindex_mat(2,i),:)=data_coils_last(i,:);
```

Pulseq – iRadon

- The beginning of the script is similar to the FID/SE/2D FFT examples
- Projection directions are calculated from the Pulseq sequence file
- Matrix size and projection data centering needs to be adjusted manually
- Example also shows a comparison between absolute value and complex inverse Radon transformations
 - Check out the gradient-echo data to notice the difference



Pulseq – 2D Gridding

- The beginning of the script is similar to the FID/SE/2D FFT examples
- Calculate the sampling trajectory from the Pulseq sequence file
 - Any crazy trajectory can be reconstructed
 - Also 2D Cartesian data are compatible
- FOW and imaging matrix size may need to be adjusted manually
- The screen shot ends with showing the gridded k-space data
- Following (not shown) part is a trivial 2D FFT and a sum-of-square combination for multichannel data

```
/home/zaitsev/pulseq_home/pulseq_public/matlab/ismrmTutorial/recon/recon04_2D_Gridding.m
[adc_len,channels,readouts]=size(data_unsorted);
rawdata = double(permute(data_unsorted, [1,3,2])); % channels last
rawdata = reshape(rawdata, [adc len*readouts,channels]);
% reconstruct the trajectory
traj recon delay=[0 0 0]*1e-6; % adjust this parameter to potentially improve resolution & geometric accuracy
                       % It can be calibrated by inverting the spiral revolution dimension and making
                       % two images match. for our Prisma and a particular trajectory we found 1.75e-6
                       % it is also possisible to provide a vector of 3 delays (varying per axis)
[ktraj_adc, t_adc, ktraj, t_ktraj, t_excitation, t_refocusing] = seq.calculateKspacePP('trajectory_delay',traj
max_abs_ktraj_adc=max(abs(ktraj_adc'));
[~, slcDim]=min(max_abs_ktraj_adc);
encDim=find([1 2 3]~=slcDim);
% figure; plot(t_ktraj, ktraj'); % plot the entire k-space trajectory
figure; plot(ktraj(encDim(1),:),ktraj(encDim(2),:),'b',...
             ktraj_adc(encDim(1),:),ktraj_adc(encDim(2),:),'r.'); % a 2D plot
axis('equal'); title('2D k-space trajectory');
Some Define FOV and resolution and simple off-resonance frequency correction
                                                                                      update for
%fov=30e-3; Nx=128; Ny=Nx; % OCRA
fov=256e-3; Nx=256; Ny=Nx; % whole-body scanners
                                                                                       your data
os=2; % oversampling factor (we oversample both in image and k-space)
offresonance=0; % global off-resonance in Hz
% rudimentary off-resonance correction
nex=length(t excitation);
if nex>1
   for e=2:nex
        il=find(t_adc>t_excitation(e),1);
            i2=max(find(t_adc<t_excitation(e+1)));
            i2=length(t adc);
        t adc ex(i1:i2)=t adc ex(i1:i2)-t excitation(e);
    rawdata(:,c) = rawdata(:,c) .* exp(-li*2*pi*t_adc_ex'*offresonance);
% here we expect Nx, Ny, deltak to be set already
% and rawdata ktraj adc loaded (and having the same dimensions)
kxm=round(os*os*Nx/2);
kym=round(os*os*Ny/2);
[kyy,kxx] = meshgrid(-kxm:(kxm-1), -kym:(kym-1));
kvv=-kvv*deltak/os: % we swap the order ind invert one sign to account for Matlab's strange column/line convent
kqd=zeros([size(kxx) channels]);
    kgd(:,:,c)=griddata(ktraj_adc(encDim(1),:),ktraj_adc(encDim(2),:),rawdata(:,c),kxx,kyy,'cubic');
kgd(isnan(kgd))=0;
figure; imagesc(log(abs(kgd(:,:,1)))); axis('square'); title('k-space data after gridding');
                                                                                                   Ln 44 Col 1
```



Outline


- Pulseq sequence definition and programming language overview
- Using *Pulseq* on different types of scanners:
 - Siemens, GE, and tabletop scanners using OCRA
- Live demo of *Pulseq* across sites
 - 1. Basic pulse-acquire scans: Free Induction Decay (FID), Spin-Echo (SE)
 - 2. Cartesian spin-eco imaging
 - 3. Radial imaging (spin-echo and gradient echo)
 - 4. Advanced topics (gradient delays and correction, etc...) are only covered in the live portion of the tutorial


How to reach us

Site	Pulseq Users & Contact Info
Vienna 7 T Siemens	Maxim Zaitsev < maxim.zaitsev@meduniwien.ac.at >
Michigan 3 T GE	Jon-Fredrik Nielsen <jfnielse@umich.edu> Melissa Haskell <mhask@umich.edu></mhask@umich.edu></jfnielse@umich.edu>
Dortmund 0.48 T	Benjamin Menkuec <benjamin.menkuec@fh-dortmund.de></benjamin.menkuec@fh-dortmund.de>
Boston 0.35 T	Lincoln Craven-Brightman: <lcraven-brightman@mgh.harvard.edu> Jason Stockmann <jstockmann@mgh.harvard.edu></jstockmann@mgh.harvard.edu></lcraven-brightman@mgh.harvard.edu>
Magdeburg 0.27 T	Marcus Prier <marcus.prier@ovgu.de> David Schote <david.schote@st.ovgu.de></david.schote@st.ovgu.de></marcus.prier@ovgu.de>
Leiden 50 mT	Tom O'Reilly <t.o_reilly@lumc.nl></t.o_reilly@lumc.nl>

Join us on Wed May 19 at 19:00 UTC

Pulseq – that's the way you do it!





...dream of a sequence in the morning...

...check the images in the afternoon!